


Elevation TPLO Plate

Locking - Polyaxial

Correct order of screw placement



Weight Plate Guidelines

NB Guidelines only.

Pre-operative planning with consideration of radiographic findings, conformation, patient factors and owner's compliance is essential.

Elevation TPLO Plates	Patient Weight
2.0mm TPLO Plate - 28mm - 6 Holes	3.5 - 7kg
2.4mm TPLO Plate - 33mm - 6 Holes	7 - 12kg
2.7mm TPLO Plate - 43mm - 6 Holes	12 - 18kg
3.5mm Small TPLO Plate - 53mm - 6 Holes	18 - 30kg
3.5mm Standard TPLO Plate - 62mm - 6 Holes	22 - 40kg
3.5mm Broad TPLO Plate - 78mm - 8 Holes	35 - 50kg
3.5mm Jumbo TPLO Plate - 88mm - 10 Holes	50kg +

Radiograph Template

Implants: Elevation Plate

Manufacturer: InVictos

